

WHEN PURITY IS PARAMOUNT

Accumax Heavy Metal Free Tubes

Accumax



Pure, Additive-free
Natural Color Caps
For Easy Labelling



Tubes And Caps
Made With Virgin,
100% Medical-grade Resins.



Certified below one part
per billion (<1ppb) trace
of heavy metals




Large writing area



Heavy Metals Free centrifuge tubes are specifically manufactured to make sure that common trace metals that can interfere with experiments are below an ICP-MS detection level of one part per billion (<1ppb).

Metals included in testing are:
Boron, Aluminum, Phosphorous,
Arsenic, Selenium, Cadmium, Lead,
Chromium, Manganese, Iron, Nickel,
Copper, Zinc.



15 ML

50 ML

For more details visit us at www.accumaximum.com