

# Plate Compatibility Chart

PLATE COMPATIBILITY CHART						
	Note/ Holder No.	APC-021	APC-022	APC-023	APC-024	APC-025
Standard Well (0.2ml, height 21mm)		✓	✓			
Low Profile Well (0.1ml, height 16mm)				✓	✓	✓
<b>Applied Biosystems</b>						
7000			✓			
7300	4367032		✓			
7500	4367032		✓			
7900HT Fast (with tray) (0.1ml type)				✓		
7900HT Fast (with tray) (0.2ml type)			✓			
StepOne	48-Well					
StepOne Plus				✓		
7500Fast	4359652			✓		
ViiA7 Fast (0.1ml model)	4459846			✓		
ViiA7 Fast (0.2ml model)	4459845		✓			
QuantStudio Fast 96-well Block (0.1ml model)	4459846			✓		
QuantStudio Fast 96-well Block (0.2ml model)	4459845		✓			
5700	Sequencer		✓			
7700	Sequencer		✓			
310/3100/3130	Sequencer		✓			
3700/3730	Sequencer		✓			
<b>Roche</b>						
LightCycler 480						✓
LightCycler 480 II						✓
LightCycler 96						✓
LightCycler Nano	32-Well					
<b>Bio-Rad</b>						
iCycler iQ		✓				
iQ5		✓				
MyiQ		✓				
Chromo4		✓				
Opticon I					✓	✓
Opticon II					✓	✓
mini Opticon	48-Well					
CFX96					✓	✓
CFX384						
<b>Stratagene</b>						
Mx4000/3000P/3005P		✓				
<b>Eppendorf</b>						
MasterCycle ep Realplex		✓				
<b>Thermo Scientific</b>						
PikoReal 24 Real-Time PCR	24-Well					

# Plate Compatibility Chart

PLATE COMPATIBILITY CHART						
	Note/ Holder No.	APC-021	APC-022	APC-023	APC-024	APC-025
<b>Biometra</b>						
Topical Thermocycler		☑				

ACCUMAX QPCR PRODUCTS CORRESPONDING TO APPLIED BIOSYSTEMS & BIO-RAD PLASTICWARE		
Accumax Cat. No.	Applied Biosystems Cat. No.	Bio-Rad Cat. No.
APC-021	N/A	MLP-9601 / 223-9441
APC-022	N8010560 / 4316813 (big pack)	MSS-9601
APC-023	4346907	N/A
APC-024	N/A	MSP-9601 / MLL-9601
APC-025	N/A	MLL-9651 (white)
APC-032	4311971	MSC-1001



**Accumax Lab Devices Pvt. Ltd.**  
www.accumaximum.com



ISO 9001:2015  
ISO 13485:2016

\* OEM request can be entertained. Specifications can be changed without notice for quality improvement.