

Plate Compatibility Chart

PLATE COMPATIBILITY CHART														
	Note/ Holder No.	APC- 018	APC- 020	APC- 021	APC- 021-TO	APC- 022	APC- 022-F	APC- 022-S	APC- 023	APC- 024	APC- 025	APC- 026	APC- 027	APC- 028
Standard Well (0.2ml, height 21mm)				✓	✓	✓	✓	✓						
Low Profile Well (0.1ml, height 16mm)			✓						✓	✓	✓			
Applied Biosystems														
7000						✓	✓	✓						
7300	4367032					✓	✓	✓						
7500	4367032					✓	✓	✓						
7900HT Fast (with tray) (0.1ml type)			✓						✓					
7900HT Fast (with tray) (0.2ml type)						✓	✓	✓						
StepOne	48-Well													
StepOne Plus			✓						✓					
7500Fast	4359652		✓						✓					
ViiA7 Fast (0.1ml model)	4459846		✓						✓					✓
ViiA7 Fast (0.2ml model)	4459845					✓	✓	✓						
QuantStudio Fast 96-well Block (0.1ml model)	4459846		✓						✓					✓
QuantStudio Fast 96-well Block (0.2ml model)	4459845					✓	✓	✓						
5700	Sequencer					✓	✓	✓						
7700	Sequencer					✓	✓	✓						
310/3100/3130	Sequencer					✓	✓	✓						
3700/3730	Sequencer					✓	✓	✓						
3730xl/3500xl												✓	✓	
Sequencing Plates												✓	✓	
3730, 3730xl, 3500, 3500xl												✓	✓	
Amersham Mega Base 4000												✓	✓	
Roche														
LightCycler 480												✓	✓	
LightCycler 480 II												✓	✓	
LightCycler 96											✓			✓
LightCycler Nano	32-Well													
Bio-Rad														
iCycler iQ				✓	✓									
iQ5				✓	✓									
MyiQ				✓	✓									
Chromo4				✓	✓									
Opticon I									✓	✓				
Opticon II									✓	✓				
mini Opticon	48-Well													
CFX96									✓	✓				✓
CFX384												✓	✓	
Stratagene														
Mx4000/3000P/3005P				✓	✓									

Note: All the plates are compatible with different End point PCR (Thermal Cycler) machines. based on respective Thermal block well number 96/384 plates can be selected and used accordingly

Plate Compatibility Chart

PLATE COMPATIBILITY CHART								
	Note/ Holder No.	APC-018	APC-021	APC-022	APC-023	APC-024	APC-025	APC-026
Eppendorf								
Master Cycle ep Realples			<input checked="" type="checkbox"/>					
Gradient Cycler ep								<input checked="" type="checkbox"/>
Thermo Scientific								
PikoReal 24 Real-Time PCR	24-Well							
Biometra								
Toptical Thermocycler			<input checked="" type="checkbox"/>					
Amersham								
Megabase 4000								<input checked="" type="checkbox"/>
Qiagen®								
Rotogene®		<input checked="" type="checkbox"/>						<input checked="" type="checkbox"/>

ACCUMAX QPCR PRODUCTS CORRESPONDING TO APPLIED BIOSYSTEMS & BIO-RAD PLASTICWARE		
Accumax Cat. No.	Applied Biosystems Cat. No.	Bio-Rad Cat. No.
APC-020	4346907	N/A
APC-021	N/A	MLP-9601 / 223-9441
APC-021-TO	N/A	MLP-9601 / 223-9441
APC-022	N8010560 / 4316813 (big pack)	MSS-9601
APC-022-F	N8010560 / 4316813 (big pack)	MSS-9601
APC-022-S	N8010560 / 4316813 (big pack)	MSS-9601
APC-023	4346907	N/A
APC-024	N/A	MSP-9601 / MLL-9601
APC-025	N/A	MLL-9651 (white)
APC-032	4311971	MSC-1001

Accumax

Accumax Lab Devices Pvt. Ltd.
www.accumaximum.com

ISO 9001:2015
ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.