

Technical Information

Power	2 x 1.2 V / 1000 mAh Rechargeable NiMH batteries.
Battery Usage Duration	Eight hours continuous use
Battery Charging Duration	2 – 3 Hours to charge fully
Volumetric Range	Plastic or Glass pipette from 0.1 ml up to 200 ml
Safety System	Replaceable hydrophobic 0.45 um PTFE filter Safety check valve
Standard Accessories	<ul style="list-style-type: none">• Battery charger• Hydrophobic PTFE membrane Filter 0.45 um – 2Pc• Hydrophobic PTFE membrane Filter 0.20 um – 1Pc• Wall mounting stand with screw• Desktop stand