

*Accumax*



IMMACULATE  
PRECISION

Accumax Reservoirs

## Accumax Reservoirs

Sophisticated procedures require manufacturing sophistication in measurement as well as storage equipments. Only this level of commitment to quality produces ground-breaking results and exceptional outcomes. Thankfully you don't have to look too far.

Accumax Reagent Reservoirs are ergonomically designed after relentless data analysis for complete ease of the user. These wide based stable reservoirs are used for multichannel pipetting. You can expect maximum retrieval and enhanced efficiency with Accumax Reservoirs.



50 ml



10 ml




25 ml

## Stability & Surety



Compact Polystyrene  
Body With Broad Edges For  
Repetitive Filling & Firm Stand




Angled Bottom  
Design For Maximum  
Reagent Recovery



Pour-off Spouts To Inhibit  
The Spillage While Pouring  
The Reagent

Graduated Wall For  
Quick Measurement Of  
The Remaining Liquid



Stable & Spill Proof  
Extra Large Base



## Ordering Information

Accumax Reservoirs		
Cat. No.	Volume	Qty/Case
ARV-10-B	10 ml Bulk	100 x 5
ARV-10-S	10 ml sterile	40 x 5
ARV-10-S-W	10 ml Individually Wrap Sterile	50
ARV-25-B	25 ml Bulk	100
ARV-25-S	25 ml sterile	40 x 5
ARV-25-S-W	25 ml Individually Wrap Sterile	100
ARV-50-B	50 ml Bulk	100
ARV-50-S	50 ml sterile	40 x 5
ARV-50-S-W	50 ml Individually Wrap Sterile	100

*Accumax*

**Accumax Lab Devices Pvt. Ltd.**  
www.accumaximum.com



ISO 9001:2015  
ISO 13485:2016

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.