

*Accumax*



A NEW DIMENSION  
TO ACCURACY

Accumax Serological Pipette

## Accumax Serological Pipette

Dispensing and transferring liquids during lab experiments were traditionally replete with errors arising due to poor graduation marks, parallax error but with time technology and innovation conquered them all. However, reaching absolute accuracy is something only a few have managed to achieve.

Thankfully the Serological Pipette by Accumax manufactured from virgin polystyrene and having a Sterility Assurance Level of  $10^{-6}$  are calibrated for heavy-duty accuracy within +/- 2% volume range. These pipettes are manufactured with care and have undergone prudent quality control to deliver the best results.



## Best-In-Class Results



Manufactured From Virgin Polystyrene In A Class 100K Clean Room




US FDA Approved Raw Material



Large & Clear Graduation Marks

Negative Graduation Marks For Extra Capacity



Certified Free From DNase, RNase, Cytotoxin & Non-pyrogenic



Calibrated For Accurate Dispensing Within +/- 2% Volume Range



## Best-In-Class Results



Filter Plug Prevents Both Fluid And Liquid Vapors From Contaminating The Pipettor And Sample



Smooth Internal Surface

Highly Transparent And Clear



Color Coding For Easy Identification



# Ordering Information

Accumax Serological Pipettes		
Cat. No.	Volume	Qty/Case
AS-1-S-W	1 ml Individually wrap Sterile	250X4
AS-2-S-W	2 ml Individually wrap Sterile	200X4
AS-5-S-W	5 ml Individually wrap Sterile	125X4
AS-10-S-W	10 ml Individually wrap Sterile	100X4
AS-25-S-W	25 ml Individually wrap Sterile	75X4
AS-50-S-W	50 ml Individually wrap Sterile	50X4
AS-100-S-W	100 ml Individually wrap Sterile	25X4

*Accumax*

**Accumax Lab Devices Pvt. Ltd.**  
www.accumaximum.com



ISO 9001:2015  
ISO 13485:2016

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.