

Accumax



FLAWLESS SIMPLICITY

Accumax Pasteur Pipette

Accumax Pasteur Pipette

Pasteur Pipettes have a long history of laboratory use and since then many advancements have been made, especially in its material and design. It becomes essential for any equipment manufacturer to improve its performance without compromising on the simplicity.

In Accumax Pasteur Pipettes you'll find that flawless design married to the simplicity of use. These non-sterile pipettes made from medical grade LDPE are excellent for drop repeatability and high precision required in a laboratory.



Qualities Crafted For Perfection



Soft And Elastic Bulb For Easy Pipetting

US FDA Approved Medical Grade Virgin Polypropylene Plastic Without Any Trace Of Metals. Highly Transparent And Clear




Manufactured In Class 100K Clean Room Facility

Smooth Orifice For Precise Dispensing



Non-toxic, Medical Grade LDPE Used

Visible External Graduation For Accurate Dispensing



Swiss-made Diamond Finish Mold For Flawless "Super Smooth" Internal Surface

Comes In Two Variants - 1 ML And 3 ML



Certified DNase, RNase, Pyrogen, Endotoxin, Human DNA & PCR Inhibitor Free.

Ordering Information

Accumax Pasteur Pipettes			
Cat. No	Volume	Quantity/Pack	Pack/Case
AP-1-B	1 ml Bulk	500	10 Packs
AP-1-B-S	1 ml Bulk Sterile	50x10 Packs	10 Packs
AP-1-S-W	1 ml Individually Wrap Sterile	500x1 Packs	10 Packs
AP-3-B	3 ml Bulk	500	10 Packs
AP-3-B-S	3 ml Bulk Sterile	50x10 Packs	10 Packs
AP-3-S-W	3 ml Individually Wrap Sterile	500x1 Packs	10 Packs

Accumax

Accumax Lab Devices Pvt. Ltd.
www.accumaximum.com



ISO 9001:2015
ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.