

*Accumax*



FLAWLESS SIMPLICITY

Accumax Transfer Pipette

## Accumax Transfer Pipette

Transfer Pipettes have a long history of laboratory use and since then many advancements have been made, especially in its material and design. It becomes essential for any equipment manufacturer to improve its performance without compromising on the simplicity.

In Accumax Transfer Pipettes you'll find that flawless design married to the simplicity of use. These non-sterile pipettes made from medical grade LDPE are excellent for drop repeatability and high precision required in a laboratory.



## Qualities Crafted For Perfection



Soft And Elastic Bulb For Easy Pipetting



Manufactured In Class 100K Clean Room Facility



Non-toxic, Medical Grade LDPE Used



Flawless mold finish with super internal smooth surface



Certified DNase, RNase, Pyrogen, Endotoxin, Human DNA & PCR Inhibitor Free.


US FDA Approved Medical Grade Virgin Polypropylene Plastic Without Any Trace Of Metals. Highly Transparent And Clear



Smooth Orifice For Precise Dispensing



Visible External Graduation For Accurate Dispensing



Comes In Two Variants - 1 ML And 3 ML



## Ordering Information

Accumax Transfer Pipette			
Cat. No	Volume	Quantity/Pack	Pack/Case
AP-1-B	1 ml Bulk	500	10
AP-1-B-S	1 ml Bulk Sterile	50x10	10
AP-1-S-W	1 ml Individually Wrap Sterile	300x1	10
AP-3-B	3 ml Bulk	500	10
AP-3-B-S	3 ml Bulk Sterile	50x10	10
AP-3-S-W	3 ml Individually Wrap Sterile	300x1	10

*Accumax*

**Accumax Lab Devices Pvt. Ltd.**  
www.accumaximum.com



ISO 9001:2015  
ISO 13485:2016

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.