

Accumax



PERFECT BLEND OF
TRANSPARENCY AND SAFETY

Accumax Cytometry Tubes

Accumax Cytometry Tubes

Made up of crystal clear polystyrene, These tubes are high quality round bottom tubes for the ease of use, reliability and versatility. These tubes are manufactured in our state-of-the-art manufacturing facility to meet rigorous quality standards. The tighter dimensional tolerances provide assurance for consistent reproducible performance in your assays.



12 X 75 MM TUBES OF 5 ML CAPACITY

Product Features




Ideal for protocols in cell isolation and flow cytometry



Dual snap cap for culturing and storage application



Withstand G - Force up to 4000g



Adaptable to wide range of flow cytometry instruments



Manufactured in human touch free class 100K clean room from US FDA approved raw material



Certified DNase, RNase, Endotoxin, Pyrogen & Human DNA free



Clear & transparent surface which doesn't hinder photometry reading



5ml Polystyrene tubes with Cap



Ordering Information

Accumax Cytometry Tubes			
Cat. No.	Description	Unit / Pack	Quantity/Case
AC-T	Flow Cytometry Tubes, Round Bottom 5 ml Tube without Cap	125 / Pack	125 X 4
AC-C	Flow Cytometry Tubes, Caps for 5ml Tube	125 / Pack	125 X 4
AC-TC	Flow Cytometry Tubes, Round Bottom 5 ml Tube with Snap Cap	125 / Pack	125 X 4
AC-TC-S-W	Flow Cytometry Tubes, Sterile - Round Bottom 5 ml Tube with Snap Cap	Individually Packed	500

Accumax

Accumax Lab Devices Pvt. Ltd.
www.accumaximum.com



ISO 9001:2015
ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.