

# PERFECT BLEND OF TRANSPARENCY AND SAFETY

*Accumax*

## Accumax Storage Tubes

Made using medical grade polystyrene (Tube) and Polyethylene (Cap)



Ideal for storage applications



Tight snap cap ensure no spillage



Withstand G - Force up to 4000g



Manufactured in human touch free class 100K clean room from US FDA approved raw material



Certified DNase, RNase, Endotoxin, Pyrogen & Human DNA free



Clear and Transparent surface



12 X 75 MM TUBES OF 5 ML CAPACITY

For more details visit us at [www.accumaximum.com](http://www.accumaximum.com)

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.