

# Plate Compatibility Chart

PLATE COMPATIBILITY CHART								
	Note/ Holder No.	APC-018	APC-021	APC-022	APC-023	APC-024	APC-025	APC-026
Standard Well (0.2ml, height 21mm)			☑	☑				
Low Profile Well (0.1ml, height 16mm)					☑	☑	☑	
<b>Applied Biosystems</b>								
7000				☑				
7300	4367032			☑				
7500	4367032			☑				
7900HT Fast (with tray) (0.1ml type)					☑			
7900HT Fast (with tray) (0.2ml type)				☑				
StepOne	48-Well							
StepOne Plus					☑			
7500Fast	4359652				☑			
ViiA7 Fast (0.1ml model)	4459846				☑			
ViiA7 Fast (0.2ml model)	4459845			☑				
QuantStudio Fast 96-well Block (0.1ml model)	4459846				☑			
QuantStudio Fast 96-well Block (0.2ml model)	4459845			☑				
5700	Sequencer			☑				
7700	Sequencer			☑				
310/3100/3130	Sequencer			☑				
3700/3730	Sequencer			☑				
3730xl/3500xl								☑
Sequencing Plates								☑
3730, 3730xl, 3500, 3500xl								☑
Amersham Mega Base 4000								☑
<b>Roche</b>								
LightCycler 480								☑
LightCycler 480 II								☑
LightCycler 96								☑
LightCycler Nano	32-Well							
<b>Bio-Rad</b>								
iCycler iQ			☑					
iQ5			☑					
MyiQ			☑					
Chromo4			☑					
Opticon I							☑	☑
Opticon II							☑	☑
mini Opticon	48-Well							
CFX96							☑	☑
CFX384								
<b>Stratagene</b>								
Mx4000/3000P/3005P			☑					

# Plate Compatibility Chart

PLATE COMPATIBILITY CHART								
	Note/ Holder No.	APC-018	APC-021	APC-022	APC-023	APC-024	APC-025	APC-026
<b>Eppendorf</b>								
Master Cycle ep Realples			<input checked="" type="checkbox"/>					
Gradient Cycler ep								<input checked="" type="checkbox"/>
<b>Thermo Scientific</b>								
PikoReal 24 Real-Time PCR	24-Well							
<b>Biometra</b>								
Toptical Thermocycler			<input checked="" type="checkbox"/>					
<b>Amersham</b>								
Megabase 4000								<input checked="" type="checkbox"/>
<b>Qiagen®</b>								
Rotogene®		<input checked="" type="checkbox"/>						<input checked="" type="checkbox"/>

ACCUMAX QPCR PRODUCTS CORRESPONDING TO APPLIED BIOSYSTEMS & BIO-RAD PLASTICWARE		
Accumax Cat. No.	Applied Biosystems Cat. No.	Bio-Rad Cat. No.
APC-021	N/A	MLP-9601 / 223-9441
APC-022	N8010560 / 4316813 (big pack)	MSS-9601
APC-023	4346907	N/A
APC-024	N/A	MSP-9601 / MLL-9601
APC-025	N/A	MLL-9651 (white)
APC-032	4311971	MSC-1001

*Accumax*

Accumax Lab Devices Pvt. Ltd.  
www.accumaximum.com



ISO 9001:2015  
ISO 13485:2016

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.