

EASY, DURABLE, ACCURATE!


Accumax

Cryovial

(Temperature Stability Upto -196°C / Volume Range 1/1.8/3.5/4.5ml)



V Bottom For Maximum Sample Recovery



Autoclavable
 121°C For 15 Mins



Better Stability With Hepta Foot For Self Standing And Storage Efficiency



Easy Sample Identification, Large And Clear Patch For Personalized Labeling




Quick Seal Closure And Leak Proof Tube Seal



Externally Threaded Without O-ring To Avoid Contamination



Certified
DNase RNase Free



For more details visit us at www.accumaximum.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.