

FILTRATION RANGE

Accumax

Spin Column Tube


2.0 ml



Receiver tubes in 2.0 ml



Certified free from RNase/DNase, human DNA, pyrogens & PCR inhibitors with sterility assurance level (SAL) of 10^{-6}



Manufactured from premium-grade autoclavable, non-wettable, virgin polypropylene



Manufactured in ISO 9001 & ISO 13485 certified 100K clean room facility



Multiple filter options like Glass Fibre filter, Polyethylene filter, Cellulose Acetate filter, Regenerated Cellulose, Nylon filter, Polyvinylidene Fluoride filter



Can Withstand up to 25000 x g fixed angle rotor and 30000 x g swing-out rotor



Includes a standard 0.8 ml round bottom filter tube



For more details visit us at www.accumaximum.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.