

*Accumax*



STURDY AND ATTESTED  
MEDIA BOTTLES

Media Bottles

PETG/PET

## Media Bottle

Material: PETG/PET

Volume: PETG : 30, 60, 125, 250, 500, 1000, 2000ml

PET : 125, 250, 500 & 1000ml

Ergonomic Media bottles are suited for biopharmaceutical processing, storage of critical reagents, Media and various buffers used in laboratories. Manufactured in a clean room facility using high-quality laboratory-grade raw materials. Packaging is designed for clean room applications. The bottle and closure are secured with a tamper-resistant heat-shrink band, the bottle body is square shaped, making it convenient for storage purposes, The surface is smooth and shiny, easy to clean, and not easy to leave a residue. Our resin materials are chosen to minimize additives and potential leachables. We don't use any fillers or plasticizers.



## Product Features




Made of USP Class VI  
Raw material, with low  
leachable and extractable

Available in various  
sizes and Material



Manufactured in  
Class 8 Clean room

Gamma Sterile  
Product




Non Pyrogenic and  
Non cytotoxic

Graduated for quick  
user visibility



\*Leakproof bottles &  
HDPE closures provided

Suitable for frozen storage  
down to -40°C to -70 °C



Heat-shrink band around closure  
& neck for tamper-evident seal

\* leak proof testing is done using water, these tests using other liquids may not yield the same results.  
Tested at 715mmhg pressure at room temp

# Ordering Information

PETG				
Cat. No.	Product Description	Pack Size	Case Size	Closure
AMB-30-PG-S	30ml Media Bottle with Cap, PETG, Sterile	24	96	HDPE
AMB-60-PG-S	60ml Media Bottle with Cap, PETG, Sterile	24	96	HDPE
AMB-125-PG-S	125ml Media Bottle with Cap, PETG, Sterile	24	48	HDPE
AMB-250-PG-S	250ml Media Bottle with Cap, PETG, Sterile	24	48	HDPE
AMB-500-PG-S	500ml Media Bottle with Cap, PETG, Sterile	12	24	HDPE
AMB-1000-PG-S	1L Media Bottle with Cap, PETG, Sterile	12	24	HDPE
AMB-2000-PG-S	2 L Media Bottle with Cap, PETG, Sterile	6	12	HDPE
AMB-30-PG-B	30ml Media Bottle with Cap, PETG, Non Sterile	40	280	HDPE
AMB-60-PG-B	60ml Media Bottle with Cap, PETG, Non Sterile	40	200	HDPE
AMB-125-PG-B	125ml Media Bottle with Cap, PETG, Non Sterile	24	96	HDPE
AMB-250-PG-B	250ml Media Bottle with Cap, PETG, Non sterile	30	60	HDPE
AMB-500-PG-B	500ml Media Bottle with Cap, PETG, Non Sterile	24	48	HDPE
AMB-1000-PG-B	1L Media Bottle with Cap, PETG, Non Sterile	12	24	HDPE
AMB-2000-PG-B	2 L Media Bottle with Cap, PETG, Non sterile	6	12	HDPE
AMB-30-PG-B-S	30ml Media Bottle with Cap, PETG, Bulk, Sterile	40	280	HDPE
AMB-60-PG-B-S	60ml Media Bottle with Cap, PETG, Bulk, Sterile	40	280	HDPE
AMB-125-PG-B-S	125ml Media Bottle with Cap, PETG, Bulk, Sterile	24	96	HDPE
AMB-250-PG-B-S	250ml Media Bottle with Cap, PETG, Bulk, Sterile	30	60	HDPE
AMB-500-PG-B-S	500ml Media Bottle with Cap, PETG, Bulk, Sterile	24	48	HDPE

PET				
Cat. No.	Product Description	Pack Size	Case Size	Closure
AMB-125-PT	125ml Media Bottle with Cap, PET, Non sterile	24	96	HDPE
AMB-250-PT	250ml Media Bottle with Cap, PET, Non sterile	30	60	HDPE
AMB-500-PT	500ml Media Bottle with Cap, PET, Non sterile	24	48	HDPE
AMB-1000-PT	1 L Media Bottle with Cap, PET, Non sterile	12	24	HDPE
AMB-125-PT-S	125ml Media Bottle with Cap, PET, Sterile	24	96	HDPE
AMB-250-PT-S	250ml Media Bottle with Cap, PET, Sterile	30	60	HDPE
AMB-500-PT-S	500ml Media Bottle with Cap, PET, Sterile	24	48	HDPE
AMB-1000-PT-S	1 L Media Bottle with Cap, PET, Sterile	12	24	HDPE

Packed in shrink-wrapped trays.

*Accumax*

**Accumax Lab Devices Pvt. Ltd.**  
www.accumaximum.com



ISO 9001:2015  
ISO 13485:2016

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.