


SAFETY CAP

ACCUMAX

Thread Size: GL 45

Available in 1, 2, 3, 4 and 5 port



High quality materials (PTFE and PP) ensure maximum safety in your laboratory



Advanced vapour valve keeps lab environment free from hazardous solvent vapour



Efficient PTFE filter for preventing contaminants mediated noise



Freely rotating cap body prevents twisting tubes during solvent bottle exchange



To be used with solvents bottles for safer lab environment and lesser uncertainty of analysis.



Chemically resistant against aggressive organic solvents



Blind fittings For closing unused ports on the cap



For more details visit us at www.accumaximum.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.