

ACCUMAX







VERSATILE, RELIABLE
AND CONSISTENT

TECAN[®] TYPE TIPS

Volume: 10, 200 and 1000 ul

TECAN® TYPE TIPS

The Accumax Robotic Tips for Tecan® Type Tips have been carefully designed, produced, and validated to ensure reliable performance when used with LiHa® Arms. These tips boast superior straightness and a low coefficient of variation, ensuring accurate and precise performance. They are compatible with the Freedom Evo® 100, Evo® 150, or Evo® 200 with Liquid Handling Arm, as well as the Freedom EVO® 75 8+1. Available in a variety of tip volumes, with or without filters, the Accumax Robotic Tips for Tecan® Type Tips are the perfect choice for your automatic liquid handling needs.



PRODUCT FEATURES




Constructed of premium-grade virgin polypropylene. Available in carbon-embedded conductive tips and clear tips




USP class VI resins



Validated on the actual instrument for which the tip is created



Polycarbonate (PC) tray with color coding to indicate different volumes



Available in Blister packs and Eco flip box

Manufactured in ISO-certified 9001 and 13485 facility




All tips are certified Endotoxin, Heavy Metals, RNase, DNase and DNA free



Nonsterile, sterile and filtered sterile formats.



Polyethylene filter to prevent aerosol contamination



Ordering Information

Conductive tips

Product Code	Description	Tray color	Per rack	Quantity/case
ATTB-10-FL	Accumax Tecan type 10µL tips Standard-Conductive, Flip box	White	2x96	2 x 96 tips x 5 Racks x 18 Packs / case
ATTB-10-F-FL	Accumax Tecan type 10µL tips Standard Filter-Conductive, Flip Box	White	2x96	2 x 96 tips x 5 Racks x 18 Packs / case
ATTB-200-FL	Accumax Tecan type 200µL tips Standard- conductive, Flip Box	Blue	2x96	2 x 96 tips x 5 Racks x 18 Packs / case
ATTB-200-F-FL	Accumax Tecan type 200µL tips Standard Filter-conductive, Flip Box	Blue	2x96	2 x 96 tips x 5 Racks x 18 Packs / case
ATTB-1000-FL	Accumax Tecan type 1000µL tips Standard-conductive, Flip box	Red	2x96	2 x 96 tips x 5 Racks x 10 Packs / case
ATTB-1000-F-FL	Accumax Tecan type 1000µL tips Standard-Filter- Conductive, Flip box	Red	2x96	2 x 96 tips x 5 Racks x 10 Packs / case
ATTB-10-BL	Accumax Tecan type 10µL tips Standard-Conductive, 1-Blister pack	White	1x96	1 x 96 tips x 24 packs / case
ATTB-10-F-BL	Accumax Tecan type 10µL tips Standard Filter-Conductive, 1-Blister pack	White	1x96	1 x 96 tips x 24 packs / case
ATTB-10-S-BL	Accumax Tecan type 10µL tips Sterile-Conductive, 1-Blister Pack	White	1x96	1 x 96 tips x 24 packs / case
ATTB-10-S-F-BL	Accumax Tecan type 10µL tips Sterile-Filter-Conductive, 1-Blister pack	White	1x96	1 x 96 tips x 24 packs / case
ATTB-200-BL	Accumax Tecan type 200µL tips Standard- conductive, 1-Blister pack	Blue	1x96	1 x 96 tips x 24 packs / case
ATTB-200-F-BL	Accumax Tecan type 200µL tips Standard Filter-conductive, 1-Blister Pack	Blue	1x96	1 x 96 tips x 24 packs / case

Product Code	Description	Tray color	Per rack	Quantity/case
ATTB-200-S-BL	Accumax tecan type 200µL tips Sterile-conductive, 1-Blister Pack	Blue	1x96	1 x 96 tips x 24 packs / case
ATTB-200-S-F-BL	Accumax Tecan type 200µL tips Sterile-Filter-Conductive, 1-Blister Pack	Blue	1x96	1 x 96 tips x 24 packs / case
ATTB-1000-BL	Accumax Tecan type 1000µL tips standard-conductive, 1-Blister pack	Red	1x96	1 x 96 tips x 24 packs / case
ATTB-1000-F-BL	Accumax Tecan type 1000µL tips Standard-Filter- Conductive, 1-Blister pack	Red	1x96	1 x 96 tips x 24 packs / case
ATTB-1000-S-BL	Accumax Tecan type 1000µL tips Sterile-conductive, 1-Blister pack	Red	1x96	1 x 96 tips x 24 packs / case
ATTB-1000-S-F-BL	Accumax Tecan type 1000µL tips Sterile-Filter-Conductive, 1-Blister pack	Red	1x96	1 x 96 tips x 24 packs / case

Clear tips

Product Code	Description	Tray color	Per rack	Quantity/case
ATTC-10-FL	Accumax Tecan type 10µL tips Standard- clear, Flip Box	White	2x96	2 x 96 tips x 5 Racks x 18 Packs / case
ATTC-10-F-FL	Accumax Tecan type 10µL tips Standard Filter- clear, Flip Box	White	2x96	2 x 96 tips x 5 Racks x 18 Packs / case
ATTC-200-FL	Accumax Tecan type 200µL tips Standard- clear, Flip Box	Blue	2x96	2 x 96 tips x 5 Racks x 18 Packs / case
ATTC-200-F-FL	Accumax Tecan type 200µL tips Standard Filter-clear, Flip Box	Blue	2x96	2 x 96 tips x 5 Racks x 18 Packs / case
ATTC-1000-FL	Accumax Tecan type 1000µL tips Standard-clear, Flip Box	Red	2x96	2 x 96 tips x 5 Racks x 10 Packs / case
ATTC-1000-F-FL	Accumax Tecan type 1000µL tips Standard-Filter-clear, Flip Box	Red	2x96	2 x 96 tips x 5 Racks x 10 Packs / case
ATTC-10-BL	Accumax Tecan type 10µL tips Standard- clear, 1-Blister Pack	White	1x96	1 x 96 tips x 24 packs / case
ATTC-10-F-BL	Accumax Tecan type 10µL tips Standard Filter- clear,1-Blister Pack	White	1x96	1 x 96 tips x 24 packs / case
ATTC-10-S-BL	Accumax Tecan type 10µL tips Sterile-clear, 1-Blister Pack	White	1x96	1 x 96 tips x 24 packs / case
ATTC-10-S-F-BL	Accumax Tecan type 10µL tips Sterile-Filter-clear, 1-Blister pack	White	1x96	1 x 96 tips x 24 packs / case
ATTC-200-BL	Accumax Tecan type 200µL tips Standard- clear,1-Blister pack	Blue	1x96	1 x 96 tips x 24 packs / case
ATTC-200-F-BL	Accumax Tecan type 200µL tips Standard Filter-clear, 1-Blister pack	Blue	1x96	1 x 96 tips x 24 packs / case

Product Code	Description	Tray color	Per rack	Quantity/case
ATTC-200-S-BL	Accumax Tecan type 200µL tips Sterile-clear, 1-Blister pack	Blue	1x96	1 x 96 tips x 24 packs / case
ATTC-200-S-F-BL	Accumax Tecan type 200µL tips Sterile-Filter-clear, 1-Blister Pack	Blue	1x96	1 x 96 tips x 24 packs / case
ATTC-1000-BL	Accumax Tecan type 1000µL tips Standard-clear, 1-Blister Pack	Red	1x96	1 x 96 tips x 24 packs / case
ATTC-1000-F-BL	Accumax Tecan type 1000µL tips Standard-Filter-clear, 1-Blister pack	Red	1x96	1 x 96 tips x 24 packs / case
ATTC-1000-S-BL	Accumax Tecan type 1000µL tips Sterile-clear, 1-Blister pack	Red	1x96	1 x 96 tips x 24 packs / case
ATTC-1000-S-F-BL	Accumax Tecan type 1000µL tips Sterile-Filter-clear, 1-Blister pack	Red	1x96	1 x 96 tips x 24 packs / case

ACCUMAX

Accumax Lab Devices Pvt. Ltd.
www.accumaximum.com

ISO 9001:2015
ISO 13485:2016

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.