

SOLUTION TO TIRED HANDS, BECAUSE PIPETTING SHOULDN'T HURT!


ACCUMAX

Accumax FAB LF Pipette




Effortless Low Force
Tip Ejection

Ergonomic design
with proper weight
distribution



The Accumax Low Force
Pipette forms a light, perfect
seal between the pipette and
tips, making tips virtually
effortless to eject.

Made using highly
durable material for
protection against
chemical corrosion &
physical shocks




Designed with a cylindrical
shaft and tip design that
provides a reliable, consistent
seal without excessive force.

UV resistant & Easy in-
house calibration &
tool free cleaning



Total hand comfort with
lower spring forces.

Highly durable
PVDF tip cone



Magnet
assisted piston for
consistently precise
results

Calibrated as per
ISO 8655-6 In ISO 17025
accredited lab



Innovative spring and
seal design

Fully
autoclavable



For more details visit us at www.accumaximum.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.